# Evaluation of the qLabs<sup>®</sup> FIB system, a novel rapid and easy-to-use point-of-care system to quantify functional fibrinogen level using a single drop of citrated whole blood sample. S. SANFILIPPO<sup>1</sup>, L. BUISSON<sup>1</sup>, H. ROUABEHI<sup>1</sup>, F. DEPASSE<sup>2</sup>, E. PEYNAUD-DEBAYLE<sup>3</sup>, B. DUMONT<sup>3</sup>, T. DONNET<sup>1</sup>

<sup>1</sup>Diagnostica Stago R&D, Gennevilliers, France; <sup>2</sup>Diagnostica Stago Clinical and Pharmaceutical Development, Asnières-sur-Seine, France; <sup>3</sup>Department of Biological Hematology and Transfusion, Louis-Mourier Hospital, Colombes, France,

### INTRODUCTION

The qLabs<sup>®</sup> FIB system includes an electrometer (Fig. 1A) and a test strip (B) for a functional fibrinogen measurement from 1.0 to 4.0 g/L range in less than 10 minutes. After citrated blood drop deposition on the strip, the meter measures blood flow rate along microfluidic thrombin-coated channels and displays fibrinogen plasma concentration adjusted by sample hematocrit calculated simultaneously.

A

### AIM

Analytical performances of the qLabs® FIB system were through evaluated а study with comparison the STA® Liquid Fib predicate assay on STA-R<sup>®</sup> Max (Stago), and a 3-sites precision study.



B&A and Passing Bablok results were compliant with specifications. Moreover, the Extreme Studentized Deviate test (ESD, Rosner 1983) did not statistically detect any outlier among the sample distribution (data not shown).

The qLabs FIB assay was demonstrated as substantially equivalent to its predicate device, the STA<sup>®</sup> Liquid Fib assay on STA-R<sup>®</sup> Max (Stago).

Electrometer reading zone Β. Timing and electrical signals collection and processing Left and right electrod pairs Double channel coated with dried thrombin Single electrod pai mixing zone 15 cm entry well Single channel coated with abciximab and polybrene. Figure 1. View of the qLabs<sup>®</sup> FIB system

Passing Bablok fit (y = 0.02 + 0.99 x)----- 95% CI 0 0,5 1 1,5 2 2,5 3 3,5 [STA-R MAX Fibrinogen] g/I Correlation Passing Bablok regression fit

Specifications		coefficient	
Specifications	Slope (a) 1 ± 0.05	Intercept 0 ± 0.10	r ≥ 0.92
Comparison method study (n=110)	0.99	0.02	0.95

Figure 4. Passing Bablok regression – qLabs<sup>®</sup> FIB vs. STA<sup>®</sup> Liquid Fib on STA-R<sup>®</sup> Max The one-site methods comparison study included 110 whole blood samples distributed over the qLabs FIB<sup>®</sup> measurement range [1 – 4 g/L]. 18% (20/110) were diluted whole blood samples (Fig. 2). For each sample, fibrinogen concentration was assessed using the qLabs<sup>®</sup> FIB and the STA<sup>®</sup> Liquid Fib assay on STA-R<sup>®</sup> Max (Stago).

A Bland & Altman and Passing Bablok distribution analyzes were applied to evaluate the substantial equivalence between qLabs FIB<sup>®</sup> and the predicate assay (EP09c directive) (Figures 3 and 4).

A 3-sites precision evaluation protocol was followed to assess the accuracy of qLabs<sup>®</sup> plasma quality controls (5 FIB using (EP05-A3 directive). The precision types were expressed with both SD and CV%, with one (repeatability/within-lab precision) calculated (reproducibility) ANOVA (Tables 1 and 2).

Site # Site # Site #

## **MATERIAL & METHODS**

#### Table 1. Within-site precision

	QC Control	Mean [Fib] (g/L)	Repeatability		Within Laboratory - Precision	
			SD	%CV	SD	%CV
1 -	Level 1	2.77	0.08	2.9%	0.08	2.9%
	Level 2	1.04	0.04	4.1%	0.04	4.1%
2 -	Level 1	2.67	0.08	2.9%	0.09	3.2%
	Level 2	0.99	0.03	2.7%	0.03	2.7%
3 -	Level 1	2.83	0.11	3.8%	0.11	3.8%
	Level 2	1.04	0.05	4.4%	0.05	4.6%

#### RESULTS

QC Control	Mean [Fib] (g/L)	Reproducibility (total 3 sites)		
		SD	%CV	
Level 1	2.76	0.12	4.4%	
Level 2	1.02	0.05	4.6%	

Repeatability/within-laboratory precision and reproducibility were satisfactory with CVs a  $\leq$  10% and  $\leq$  8% respectively, according to acceptance criteria.

# CONCLUSION

The qLabs FIB<sup>®</sup> system provides a reliable and precise measurement of fibrinogen concentration from a single drop of citrated whole blood.

Further clinical trials should confirm its ability to guide management of critical bleeding, especially in the evaluation of postpartum hemorrhage severity and eligibility for fibrinogen infusion.



#### ucibility

