

st|hemE
by Stago

eSolutions
are in our blood.



eSolutions
CATALOGUE


Stago

Diagnostics is in our blood.

This document is intended for informational purposes only.


This Catalogue has no legal status of its own - it is an informative document and should be used as such.




SUMMARY

eSolutions - sthemE

stheme Connect	4
stheme Manager	5
stheme Expert QC	6
stheme Qualiris	7





sthemE Connect

sthemE Connect serves as the cornerstone of Stago digital services bundle: automatically and securely transmit your sthemE Expert QC¹ data and CPRR² reports to Stago servers. Receive optimized troubleshooting through remote support sessions and monitoring of your Stago solutions.

To order³

Please consult Stago or one of its official distributors.



¹ This service is already available under the name My Expert QC


² Cost Per Reportable Result

³ Availability depending upon the country

sthemE

by Stago

Manager



sthemE Manager

sthemE Manager is the Stago Data Manager, your cockpit to pilot your lab daily activity: result auto-validation, TAT alerts, Stago patented expert rules, full Quality Control management including Westgard rules, method validation and centralized maintenance dashboard.

To order*

Please consult Stago or one of its official distributors.

86661 sthemE Manager Starter Pack Max

86662 sthemE Manager Starter Pack sthemO

86665 Workcell license

86675 Additional site license

86676 Additional Max connection license

86677 Additional sthemO connection license

86666 sthemE Manager server license - pre-production

89465 sthemE Connect Gateway kit

The screenshot displays the 'sthemE Manager' interface with a 'PATIENTS > WORK LIST' header. It features a summary bar with statistics: 17 Order received, 5 In progress, 24 Run, 9 To be validated, 0 Validated, and 3 Transmitted. Below this is a table of patient data with columns for Sample ID, Status, Origin, Name, First name, Demand date, and Back to Position. The table lists various patients with their respective sample IDs and status.

Sample ID	Status	Origin	Name	First name	Demand date	Back to Position	Center
1007210603_03	Run	001 - ASN	AMERICAN	Christopher	05/06/2021 08:29	200994/1	COAB
1007210610_07	To be validated	001 - ASN	AUSTRALIAN	Gabriel	18/06/2021 08:09	0001 - 0001_2	200994/1
1007210612_05	Order received	001 - ASN	AUSTRALIAN	Concor	13/06/2021 08:03	200994/1	COAB
1007210612_23	Order received	001 - ASN	AMERICAN	Margaret	12/06/2021 08:05	200994/1	COAB
1007210613_01	Run	001 - ASN	AUSTRALIAN	Alexandra	26/06/2021 11:21	0001 - 0001_3	200994/1
1007210613_02	Run	001 - ASN	ENGLAND	Lily	26/06/2021 11:21	0001 - 0001_2	200994/1
1007210613_03	Run	001 - ASN	AMERICAN	Christopher	26/06/2021 11:21	0001 - 0001_3	200994/1
1007210613_05	Run	001 - ASN	AUSTRALIAN	Concor	26/06/2021 11:22	0001 - 0001_2	200994/1
1007210613_06	Order received	001 - ASN	ENGLAND	Madeline	26/06/2021 11:22	200994/1	COAB
1007210613_07	To be validated	001 - ASN	AUSTRALIAN	Gabriel	26/06/2021 11:22	0001 - 0001_3	200994/1
1007210620_01	Order received	001 - ASN	DOXALL	Alexandra	26/06/2021 14:06	200994/1	COAB
1007210620_02	Order received	001 - ASN	DOBSON	Lily	26/06/2021 14:06	200994/1	COAB
1007210620_03	Order received	001 - ASN	EHLERS	Christopher	26/06/2021 14:06	200994/1	COAB
1007210620_04	Order received	001 - ASN	GARDNER	Chelsea	26/06/2021 14:06	200994/1	COAB

* Availability depending upon the country

sthemE Expert QC

sthemE Expert QC¹ is a robust IQC peer-comparison tool: effortlessly assess the quality of your entire coagulation menu, using an automatic upload process through sthemE Manager. sthemE Expert QC helps identify potential abnormal conditions before they occur.

To order²

Please consult Stago or one of its official distributors.

82127 sthemE Expert QC 1-year contract



¹ This service is also available under the name My Expert QC

² Availability depending upon the country

sthemE Qualiris

Combined with your internal quality control, an External Quality Assessment (EQA) program will allow a comparison of performances with other participants.

Unassayed plasma covering several concentration levels and reflecting most clinical settings are provided for a regular testing. Peer group comparisons with color-coded alerts and graphs are provided few days after the data submission deadline in monthly reports. Annual reports allow month to month comparison of the results.

sthemE Qualiris can be used by all laboratories, regardless of reagents or instrumentation used.


To order*

Please consult Stago or one of its official distributors.

- 01044** Qualiris QC Premium S1
- 01045** Qualiris QC Premium S2
- 01046** Qualiris QC Experience
- 85093** Qualiris QC Premium Factors
- 85094** Qualiris QC Premium Thrombophilia
- 01047** Qualiris QC Heparin HBPM/LMWH
- 01048** Qualiris QC Heparin HNF/UFH
- 01049** Qualiris QC D-Dimer
- 01063** Qualiris QC Lupus Anticoagulant
- 01272** Qualiris QC DOAC Anti-Xa
- 01187** Qualiris QC Dabigatran
- 01188** Qualiris QC FM-FDP
- 01536** Qualiris HIT Activity
- KT0048** Qualiris Quantra QPlus
- 01050** Qualiris Diagnostic Challenge

*Availability depending upon the country

YOUDEN PLOT FOR ONE PARAMETER MONTHLY REPORT



sthemE is a trademark of the Stago Group.
The rights of the trademarks and logos used in this document belong to the Stago Group.
The use of these trademarks is not permitted without permission from the Stago Group.



sthemE
by Stago

Diagnostica Stago S.A.S.
3 allée Thérèse - CS 10009
92600 Asnières sur Seine Cedex - France
+33 (0)1 46 88 20 20
+33 (0)1 47 91 08 91
webmaster@stago.com
www.stago.com